

# Peakjuggler Manual

Doblmann Johannes

May 23, 2016

*Peakjuggler* is a node for Proteome Discoverer that aims to compare areas of peptides over several runs. Additionally it tries to find peptides that aren't identified in all runs.

## 1 Installation

### 1.1 Dependencies

- **Proteome Discoverer** version 2.0 or 2.1 is required
- **R** For some parts an R installation is needed. It is included in the installer

### 1.2 Installation

- Download the installer from <http://ms.imp.ac.at?goto=peakjuggler>

Peakjuggler consists of two nodes, one for the processing workflow and one for the consensus workflow.

## 2 Usage

### 2.1 Parameters

#### 2.1.1 Processing Node

##### 1. Search Parameters

- **Confidence Level:** Only use PSMs with at least this confidence
- **Minimum Score:** The minimum score a PSM needs to be used. This is dependent on the search engine you use

| Parameters                       |                                    |
|----------------------------------|------------------------------------|
| Hide Advanced Parameters         |                                    |
| ▼ 1. Search Parameters           |                                    |
| Confidence Level                 | Medium                             |
| Minimum Score                    | 100                                |
| Minimum Sequence Length          | 7                                  |
| Search Engine Rank               | 1                                  |
| Mass Tolerance                   | 5 ppm                              |
| ▼ 2. PhosphoRS/ptmRS Settings    |                                    |
| Probability Threshold            | 75                                 |
| PhosphoRS Column Name            | ptmRS [4]: Best Site Probabilities |
| ▼ 3. RT Correction Parameters    |                                    |
| RT Correction                    | True                               |
| ▼ 4. Performance Parameters      |                                    |
| Workpackage size                 | 10000                              |
| ▼ 5. Area Calculation Parameters |                                    |
| Minimum Width for Peak (in min)  | 0.05                               |
| Noise Level                      | 10000                              |

(a) Processing Workflow

| Parameters                        |                  |
|-----------------------------------|------------------|
| Show Advanced Parameters          |                  |
| ▲ 1. Protein Area                 |                  |
| Usage of peptides                 | average          |
| Peptides to use for Protein Area: | top x per sample |
| Use shared peptides for           | protein groups   |
| # peps for protein aera           | 3                |
| ▲ 2. Peptide Area                 |                  |
| Minimum Area                      | 50000            |

(b) Consensus Workflow

Figure 1: Node Parameters

- **Minimum Sequence Length:** The minimum number of AAs a peptide needs to have to be used
  - **Search Engine Rank:** Up to which search engine rank PSMs are considered
2. PhosphoRS/ptmRS Settings
    - **Probability Threshold:** Everything above this threshold is counted as high
    - **PhosphoRS Column:** The name of the phosphoRS/ptmRS column. If this name is not found, PJ tries some standard names, if still no column is found, modifications are taken from the search engine
  3. RT Correction Parameters
    - **RT Correction:** Peakjuggler can correct retention time shifts. This parameter activates/deactivates this feature.
  4. Performance Parameters
    - **Workpackage size:** This regulates the number of spectra that are read into the RAM. If you experience that your RAM is filled up during analysis, it is a good idea to set this parameter lower. Unfortunately with lower package size the analysis will take longer.
  5. Area Calculation Parameters
    - **Minimum Width for Peaks:** The minimum time for a peak
    - **Noise Level:** All signals below this intensity are considered as zero
    - **Mass Tolerance:** Allowed mass deviation in the XIC

## 2.1.2 Consensus Node

1. Protein Area

- **Usage of peptides:** Which method should be used to combine peptide areas to the protein area. The possibilities are:
  - sum
  - average
  - median
- **Peptides to use for Protein Area:**
  - Top x per sample: Takes the top x peptides per sample and performs the chosen method on their areas.
  - Top x over all samples: This method first sums all areas of each peptide in all samples and then takes the top x to calculate the protein area. This ensures that the same peptides are selected in all samples.
- **Use shared peptides for:**
  - All proteins: Shared peptides contribute to all proteins they appear in
  - Nothing: Shared peptides are ignored
  - Protein Groups: If a peptide is only shared within the same protein group use it
- **# peps for protein area:** The maximum number of peptides to combine for the protein area. If a protein has less peptides than this number, all available are used.

## 2. Peptide Area

- **Minimum Area:** Peptides with a smaller area are counted as not quantified

## 2.2 Workflow

### 2.2.1 Processing Node

The processing node needs two connections, one to the *Spectrum Files* node, and another one to the PSM Validation node (e.g. *Percolator* or *Target Decoy PSM Validator*).

### 2.2.2 Consensus Node


The consensus node needs to be connected to the *Protein Grouping* node.

## 3 Results


### 3.1 New columns

#### 3.1.1 Protein tab

The protein table receives two new columns: *Peakjuggler Area* and *Identified by*. The areas are colour coded to give a quick overview if the areas are similar in the samples. The



(a) Processing Workflow



(b) Consensus Workflow

Figure 2: Sample Workflows

*Identified by* column has one box per sample and indicates the origin of quantification. (MS/MS, MBR or none) See figure 3

| Proteins |                          | Protein Groups         | Peptide Groups | PSMs              | MS/MS Spectrum Info | PJ Features              | QuanResult |  |
|----------|--------------------------|------------------------|----------------|-------------------|---------------------|--------------------------|------------|--|
| Year     | Checked                  | Protein FDR Confidence | Master         | Peakjuggler Areas |                     | Identified by            | Accession  | Description  |
| 2013     | <input type="checkbox"/> | ●                      | ✓              | 2.1e8             | 2.2e8               | <input type="checkbox"/> | 51338637   | RecName: Full=60S ribosomal protein L23a                   |
| 2014     | <input type="checkbox"/> | ●                      | ✓              | 2.0e7             | 2.2e7               | <input type="checkbox"/> | 2495344    | RecName: Full=Heat shock protein 105 kDa; AltName: Full=   |
| 2015     | <input type="checkbox"/> | ●                      | ✓              | 1.4e6             | 1.4e6               | <input type="checkbox"/> | 55584151   | RecName: Full=Glucosylceramidase; AltName: Full=Acid be    |
| 2016     | <input type="checkbox"/> | ●                      | ✓              | 1.5e7             | 1.5e7               | <input type="checkbox"/> | 296439328  | RecName: Full=Importin subunit alpha-5; AltName: Full=Ka   |
| 2017     | <input type="checkbox"/> | ●                      | ✓              | 6.1e5             | 6.1e5               | <input type="checkbox"/> | 29840826   | RecName: Full=Nardilysin; AltName: Full=N-arginine dibasi  |
| 2018     | <input type="checkbox"/> | ●                      | ✓              | 1.6e7             | 1.9e7               | <input type="checkbox"/> | 74743691   | RecName: Full=Heterochromatin protein 1-binding protein 3  |
| 2019     | <input type="checkbox"/> | ●                      | ✓              | 1.5e8             | 1.5e8               | <input type="checkbox"/> | 74749006   | RecName: Full=Putative 60S ribosomal protein L13a protein  |
| 2020     | <input type="checkbox"/> | ●                      | ✓              | 3.6e7             | 3.7e7               | <input type="checkbox"/> | 3023628    | RecName: Full=ATP-dependent RNA helicase DDX3X; AltN       |
| 2021     | <input type="checkbox"/> | ●                      | ✓              | 6.7e5             | 6.4e5               | <input type="checkbox"/> | 74758580   | RecName: Full=WD repeat-containing protein 82; AltName:    |
| 2022     | <input type="checkbox"/> | ●                      | ✓              | 4.6e6             | 4.7e6               | <input type="checkbox"/> | 13124471   | RecName: Full=Ras-related protein Rab-9B; AltName: Full=   |
| 2023     | <input type="checkbox"/> | ●                      | ✓              | 2.1e6             | 2.2e6               | <input type="checkbox"/> | 27805669   | RecName: Full=Pantothenate kinase 4; Short=hPank4; AltN    |
| 2024     | <input type="checkbox"/> | ●                      | ✓              | 1.9e6             | 2.1e6               | <input type="checkbox"/> | 13124451   | RecName: Full=Ribonucleases P/MRP protein subunit POP      |
| 2025     | <input type="checkbox"/> | ●                      | ✓              | 9.7e5             | 9.8e5               | <input type="checkbox"/> | 23503074   | RecName: Full=Thymidine kinase, cytosolic                  |
| 2026     | <input type="checkbox"/> | ●                      | ✓              | 3.6e5             | 4.9e5               | <input type="checkbox"/> | 1345650    | RecName: Full=Collagen alpha-1(VII) chain; AltName: Full=  |
| 2027     | <input type="checkbox"/> | ●                      | ✓              | 6.2e6             | 7.3e6               | <input type="checkbox"/> | 74750895   | RecName: Full=Leucine-rich repeat-containing protein 47    |
| 2028     | <input type="checkbox"/> | ●                      | ✓              | 8.9e6             | 9.3e6               | <input type="checkbox"/> | 117949389  | RecName: Full=CD109 antigen; AltName: Full=150 kDa TG      |
| 2029     | <input type="checkbox"/> | ●                      | ✓              | 3.5e7             | 3.2e7               | <input type="checkbox"/> | 76803555   | RecName: Full=Nucleolar RNA helicase 2; AltName: Full=D    |
| 2030     | <input type="checkbox"/> | ●                      | ✓              | 1.8e7             | 1.6e7               | <input type="checkbox"/> | 24638339   | RecName: Full=Vesicle-associated membrane protein-asso     |
| 2031     | <input type="checkbox"/> | ●                      | ✓              | 4.5e6             | 5.1e6               | <input type="checkbox"/> | 47606219   | RecName: Full=Multifunctional methyltransferase subunit TI |
| 2032     | <input type="checkbox"/> | ●                      | ✓              | 1.9e7             | 1.9e7               | <input type="checkbox"/> | 74735389   | RecName: Full=U4/U6.U5 tri-snRNP-associated protein 1; L   |
| 2033     | <input type="checkbox"/> | ●                      | ✓              | 1.5e6             | 1.9e6               | <input type="checkbox"/> | 67476671   | RecName: Full=Oxygen-dependent coproporphyrinogen-III u    |
| 2034     | <input type="checkbox"/> | ●                      | ✓              | 3.0e7             | 3.3e7               | <input type="checkbox"/> | 2498464    | RecName: Full=28 kDa heat- and acid-stable phosphoprote    |
| 2035     | <input type="checkbox"/> | ●                      | ✓              | 5.2e7             | 5.4e7               | <input type="checkbox"/> | 5921197    | RecName: Full=Calumenin; AltName: Full=Crocalbin; AltNa    |
| 2036     | <input type="checkbox"/> | ●                      | ✓              | 5.9e6             | 6.0e6               | <input type="checkbox"/> | 313104206  | RecName: Full=Palladin; AltName: Full=SIH002; AltName: f   |
| 2037     | <input type="checkbox"/> | ●                      | ✓              | 4.5e7             | 4.5e7               | <input type="checkbox"/> | 38258929   | RecName: Full=DNA-dependent protein kinase catalytic su    |
| 2038     | <input type="checkbox"/> | ●                      | ✓              | 2.1e8             | 2.1e8               | <input type="checkbox"/> | 125969     | RecName: Full=40S ribosomal protein S4; AltName: Full=3;   |
| 2039     | <input type="checkbox"/> | ●                      | ✓              | 2.1e7             | 2.1e7               | <input type="checkbox"/> | 47117647   | RecName: Full=Actin-related protein 3; AltName: Full=Actin |
| 2040     | <input type="checkbox"/> | ●                      | ✓              | 8.0e5             | 7.4e5               | <input type="checkbox"/> | 8928568    | RecName: Full=Tumor suppressor p53-binding protein 1; SI   |
| 2041     | <input type="checkbox"/> | ●                      | ✓              | 4.8e5             | 5.2e5               | <input type="checkbox"/> | 269849532  | RecName: Full=Aldehyde dehydrogenase family 16 membe       |
| 2042     | <input type="checkbox"/> | ●                      | ✓              | 1.2e7             | 1.2e7               | <input type="checkbox"/> | 223634676  | RecName: Full=Protein TFG; AltName: Full=TRK-fused gen     |

Figure 3: Protein Table

### 3.1.2 Peptide Groups tab

Exactly the same as in the protein table. See figure 4

## 3.2 New tables

### 3.2.1 PjFeatures

This table is similar to the *PSMs* table, but not on PSM but on PCM (peptide charge modification) level (see fig 5).

It also provides the user with information regarding the quantification like peak start and end times and a button to show the extracted ion chromatogram (XIC) including the peak boundaries (see fig. 6)

### 3.2.2 PjQuanResult

The *QuanResult* is the equivalent to the *Peptide Groups* table for Features. (see fig. 7) It also features an XIC view over all samples (see fig. 8)

| Proteins | Protein Groups           | Peptide Groups                      | PSMs          | MS/MS Spectrum Info | PJ Features                    | QuanResult  |          |
|----------|--------------------------|-------------------------------------|---------------|---------------------|--------------------------------|---|----------|
| Checked  | Confidence               | PeakJugger Areas                    | Identified by | Sequence            | Modifications                  | Quality PEP   |          |
| 28       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.1e6         | 2.8e6               | DPAAQMSQGEATQSGAR              |   | 0.000136 |
| 29       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 6.0e6         | 9.7e6               | PVCAIQMLSEMIK                  | 1xCarbamidomethyl [C4]; 2xOxidation [M8(100); M12(100)] | 6.22e-07 |
| 30       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.5e7         | 2.6e7               | VESLDVDSSEAK                   |   | 7.02e-05 |
| 31       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 5.1e6         | 6.2e6               | IPPELLLRPR                     |   | 3.72e-05 |
| 32       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.2e7         | 1.7e7               | VCNFLASQVPPPSR                 | 1xCarbamidomethyl [C2]                                  | 5.68e-06 |
| 33       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 6.9e5         | 9.3e5               | GHTVTEPIQLEPELPGEGQPEAR        |   | 8.35e-06 |
| 34       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.1e7         | 1.2e7               | HRDYETATLSDIK                  |   | 1.19e-06 |
| 35       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.3e7         | 1.4e7               | LQTLVSEQPNK                    |   | 8.05e-05 |
| 36       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 7.6e6         | 9.2e6               | DGTILCTMNNK                    | 1xCarbamidomethyl [C6]                                  | 3.72e-05 |
| 37       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 4.3e6         | 5.4e6               | SYDFEFMVEK                     |   | 4.08e-06 |
| 38       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 5.2e6         | 4.2e6               | AVGASFLYEPK                    |   | 0.000115 |
| 39       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.8e7         | 3.2e7               | DYAFVHFEDR                     |   | 5.16e-05 |
| 40       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 4.8e6         | 3.9e6               | SAYQEAMDISK                    | 1xOxidation [M7(100)]                                   | 2.08e-05 |
| 41       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.7e7         | 1.9e7               | FNECGHVLADIK                   | 1xCarbamidomethyl [C4]                                  | 4.66e-07 |
| 42       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 3.8e6         | 4.1e6               | LVTQDTENELK                    |   | 6.91e-06 |
| 43       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.9e7         | 2.0e7               | KEELGAALAK                     |   | 4.7e-06  |
| 44       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.4e7         | 1.4e7               | FAEALGSTEAK                    |   | 1.66e-05 |
| 45       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 4.2e6         | 4.1e6               | KNPEVPVNFVFAESK                |   | 2.6e-07  |
| 46       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.7e6         | 3.2e6               | QLAEQEELER                     |   | 6.87e-05 |
| 47       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 5.2e6         | 7.6e6               | GLLGCLNPLQR                    | 1xCarbamidomethyl [C5]                                  | 2.21e-05 |
| 48       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.5e6         | 1.7e6               | LVQAPLDADDGNVLQEK              |   | 1.94e-06 |
| 49       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.2e6         | 2.2e6               | VVDALGNADIGK                   |   | 6.33e-06 |
| 50       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.1e6         | 2.0e6               | HQAFEALSNAGSR                  |   | 3.86e-07 |
| 51       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.9e6         | 2.4e6               | NQGFDDVLDVDTAGR                |   | 2.07e-06 |
| 52       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 6.4e6         | 6.0e6               | AVGEIVQDLYDSDK                 |   | 7.82e-06 |
| 53       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.4e7         | 1.4e7               | GTYLATFHQR                     |   | 6.33e-06 |
| 54       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.3e8         | 2.3e8               | YGVSGYPTLK                     |   | 5.47e-05 |
| 55       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 4.7e6         | 9.8e6               | LLAAVAATPPDAPNREEVFDER         |   | 7.45e-06 |
| 56       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.1e6         | 7.3e5               | TSPADHGGVSGSEGGSAVDSVAGEHSVSGR |   | 5.93e-06 |
| 57       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.4e7         | 1.7e7               | QADLYSEGLHPR                   |   | 2.48e-06 |
| 58       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 3.0e7         | 3.9e7               | VTEYGLTFTEK                    |   | 9.2e-06  |
| 59       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 2.6e7         | 2.6e7               | LAILGHNEVSK                    |   | 1.25e-06 |
| 60       | <input type="checkbox"/> | <input checked="" type="checkbox"/> | 1.0e7         | 1.1e7               | DADPLLISLR                     |   | 5.96e-05 |

Figure 4: Peptide Groups table

| Proteins                 | Protein Groups    | Peptide Groups | PSMs      | MS/MS Spectrum Info | PJ Features | QuanResult |  |          |  |             |               |        |                 |        |            |               |
|--------------------------|-------------------|----------------|-----------|---------------------|-------------|------------|--|----------|--|-------------|---------------|--------|-----------------|--------|------------|---------------|
| Checked                  | Sequence          | Peak Start     | Peak Apex | Peak End            | Charge      | MZ         | Modifications                                  | Area     | Spectrum File                              | MasterScan% | Identified By | Reason | TheoreticalMass | Score  | Confidence | Show Spectrum |
| <input type="checkbox"/> | AKQVDVDEYFFAK     | 104.80         | 104.36    | 104.62              | 2           | 620.3684   |  | 1.88E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1638.77510      | +0.262 | High       | Show Spectrum |
| <input type="checkbox"/> | ESLEALLQR         | 118.75         | 119.02    | 119.30              | 2           | 629.7920   |  | 9.22E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1059.68605      | -0.059 | High       | Show Spectrum |
| <input type="checkbox"/> | RPELLTHSTVEYQRR   | 56.43          | 56.73     | 57.08               | 4           | 467.00290  |  | 7.07E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1864.98795      | -0.297 | High       | Show Spectrum |
| <input type="checkbox"/> | RPELLHSTHVEYQRR   | 57.06          | 57.39     | 57.73               | 4           | 467.00290  |  | 7.28E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MBR           | 0      | 1864.98795      | -0.302 | High       | Show Spectrum |
| <input type="checkbox"/> | NANVWMEYRK        | 56.56          | 56.72     | 56.91               | 2           | 586.70394  |  | 2.19E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1163.53030      | -0.183 | High       | Show Spectrum |
| <input type="checkbox"/> | ESLEALLQR         | 118.97         | 119.21    | 119.41              | 2           | 629.79437  |  | 4.25E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1059.68605      | -0.088 | High       | Show Spectrum |
| <input type="checkbox"/> | AVDALPPPCESACTVPI | 118.61         | 119.27    | 120.21              | 3           | 785.71075  | C10(Carbamidomethyl)High, C14(Carbamidomethyl) | 1.66E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 2385.12067      | -0.170 | High       | Show Spectrum |
| <input type="checkbox"/> | AVDALPPPCESACTVPI | 118.84         | 119.44    | 120.25              | 3           | 785.71082  | C10(Carbamidomethyl)High, C14(Carbamidomethyl) | 1.39E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 2385.12067      | -0.252 | High       | Show Spectrum |
| <input type="checkbox"/> | RQGEAPLPIK        | 55.53          | 55.58     | 57.33               | 2           | 519.29541  |  | 1.52E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1035.50332      | -0.239 | High       | Show Spectrum |
| <input type="checkbox"/> | RQGEAPLPIK        | 55.93          | 56.37     | 56.70               | 2           | 519.29541  |  | 1.48E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1035.50332      | -0.008 | High       | Show Spectrum |
| <input type="checkbox"/> | VTVOTGVPAASEEK    | 82.29          | 82.60     | 82.87               | 2           | 722.87854  |  | 1.57E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1444.75297      | -0.300 | High       | Show Spectrum |
| <input type="checkbox"/> | AVGPGQLHGGVGR     | 82.09          | 82.38     | 82.78               | 3           | 426.23444  |  | 3.29E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1276.69991      | -0.389 | High       | Show Spectrum |
| <input type="checkbox"/> | AVGPGQLHGGVGR     | 82.38          | 82.65     | 83.10               | 3           | 426.23502  |  | 3.13E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1276.69991      | -0.186 | High       | Show Spectrum |
| <input type="checkbox"/> | EEDSLSLDGADSTGVV  | 99.10          | 99.37     | 99.63               | 2           | 925.43469  |  | 1.11E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1649.86516      | -0.163 | High       | Show Spectrum |
| <input type="checkbox"/> | YNEETFGYEVPIK     | 118.80         | 119.07    | 119.49              | 2           | 794.88269  |  | 9.70E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1568.75297      | -0.102 | High       | Show Spectrum |
| <input type="checkbox"/> | YNEETFGYEVPIK     | 118.93         | 119.26    | 119.76              | 2           | 794.88098  |  | 8.81E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1568.75297      | -0.343 | High       | Show Spectrum |
| <input type="checkbox"/> | VTVOTGVPAASEEK    | 81.87          | 82.35     | 82.63               | 2           | 722.88019  |  | 1.99E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1444.75297      | -0.243 | High       | Show Spectrum |
| <input type="checkbox"/> | HWMLDK            | 56.90          | 57.62     | 58.77               | 2           | 415.20547  |  | 4.10E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MBR           | 0      | 829.40252       | -0.294 | High       | Show Spectrum |
| <input type="checkbox"/> | HWMLDK            | 56.26          | 56.99     | 58.21               | 2           | 415.20547  |  | 4.47E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 829.40252       | -0.359 | High       | Show Spectrum |
| <input type="checkbox"/> | GELLALVK          | 102.95         | 102.78    | 102.93              | 2           | 421.77176  |  | 3.20E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 842.53458       | -0.192 | High       | Show Spectrum |
| <input type="checkbox"/> | PSHRDQVEEGEDDK    | 99.84          | 99.10     | 99.36               | 2           | 967.35773  |  | 1.99E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1723.71368      | -0.056 | High       | Show Spectrum |
| <input type="checkbox"/> | LDLFANVYHK        | 118.70         | 119.03    | 119.22              | 2           | 627.86383  |  | 3.27E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MS/MS         | 0      | 1254.72048      | -0.219 | High       | Show Spectrum |
| <input type="checkbox"/> | LDLFANVYHK        | 119.04         | 119.25    | 119.45              | 2           | 627.86383  |  | 2.55E+06 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1254.72048      | -0.067 | High       | Show Spectrum |
| <input type="checkbox"/> | TVSOLDQK          | 71.75          | 72.21     | 72.91               | 2           | 509.77448  |  | 2.24E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_02.raw | 0           | MBR           | 0      | 1018.54152      | -0.138 | High       | Show Spectrum |
| <input type="checkbox"/> | EVDEGMVIVQVK      | 55.57          | 57.30     | 56.75               | 2           | 731.94883  | 1M(Oxidation)High                              | 4.29E+07 | 20191127_QexHF1_RSLCS_helo01ug_1711_01.raw | 0           | MS/MS         | 0      | 1462.69423      | -0.362 | High       | Show Spectrum |

Figure 5: Feature table






Figure 8: XIC of a Feature in different samples

### 3.3 Scoring

Peakjuggler also tries to give the quantification a score which indicates how well the integration went. Base for this score was the DeMix-Q paper by Zhang et al.<sup>1</sup> Peakjuggler also performs a target-decoy quantification and assigns confidences to the Features. **This feature is still under development!**

<sup>1</sup>Bo Zhang, Lukas Kall, and Roman A. Zubarev

DeMix-Q: Quantification-centered Data Processing Workflow

Mol Cell Proteomics mcp.O115.055475. First Published on January 4, 2016,

doi:10.1074/mcp.O115.055475